

ELISpot PVDF plates, white

Product code:	3654-WP-10
Contents:	Contains 10 white, sterile PVDF membrane ELISpot plates (MSIPS4W10 from Millipore).
Applications:	ELISpot assay
Instructions:	To obtain maximal antibody binding capacity, the plates need to first be pre-treated, i.e. activated, with ethanol. Remove the ELISpot plate from the package and activate the membrane by adding 35% ethanol, 15 μ l per well for maximum 1 minute. Wash plate 5 times with sterile water, 200 μ l/well, and proceed according to ELISpot protocol. Once activated, it is essential that the membrane is not allowed to dry during the ELISpot procedure.
Storage:	Shipped at ambient temperature. On arrival store plates at room temperature.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Distribuito in ITALIA da
Li StarFish S.r.l.
Via Cavour, 35
20063 Cernusco S/N (MI)
telefono 02-92150794
fax 02-92157285
info@listarfish.it
www.listarfish.it