

For research use only
Cat. No. IP11072 | 48 Tests

Chicken Anemia Virus Detection Kit

Test for the detection of Chicken Anemia Virus by one-step PCR

User Manual

REV.2.2

iNtRON Biotechnology



iNtRON Biotechnology, Inc.

Phone: +82-505-550-5600 | Fax: +82-505-550-5660

e-mail: info@intronbio.com | www.intronbio.com

Table of Contents

1.	Description	1
2.	Storage	1
3.	Contents	1
4.	Specimen	1
5.	Additional required materials	2
6.	Procedure	2
6.1	DNA preparation	2
6.2	Amplification	2
6.3	Detection of amplification product	3
6.4	Interpretation	3
6.5	Elimination of carry-over contamination	3
7.	Notice	4
8.	Trouble shooting	4
9.	Ordering information	4

1. DESCRIPTION

The chicken anemia virus (CAV) virus is a 25 nm, nonenveloped, icosahedral virus with a single-stranded, circular DNA genome. The proposed classification of CAV is a new family of viruses, designated Circoviridae. The disease has been described in most countries with a developed chicken industry.

The primary importance of CAV in commercially produced chickens is due to the immunosuppression it causes in young chickens. Immunosuppression can adversely affect production characteristics including decreased weight gain, increased unevenness, poor response to vaccination and increased condemnations due to secondary infection.

CAV confirmation requires detection of virus antigens or DNA in the thymus or bone marrow. Virus isolation also can be used but is slow and expensive.

Chicken anemia virus Detection Kit is direct detection of chicken anemia virus on the basis of a genetic database, so it can diagnose very fast and accurately. It can amplify only specific gene using the PCR (Polymerase Chain Reaction) method, and take only 2~3 hours for detection. Therefore, it is a very fast accurate, reliable technique.

2. STORAGE

The components of **Chicken anemia virus Detection Kit** should be stored at -20℃, under this condition, the kit is stable until expiration date stated on the label.

3. CONTENTS

Chicken anemia virus PCR Pre-mixture	48 tubes
DNase/RNase-free water (white cap)	1 vial
CAV positive control (Yellow cap)	2 vial

Component in 20μl reaction

i-StarTaq™ DNA Polymerase
dNTPs
PCR Reaction buffer
Chemical stabilizer
Gel loading buffer
8-MOP (dissolved in DMSO)
Primers for CAV

4. SPECIMEN

Performs the test with whole blood, thymus, bone marrow, spleen, liver or feces. The specimen should be stored at -20℃ prior to use.

5. ADDITIONAL REQUIRED MATERIALS

- Disposable gloves
- DNA extraction kit (see 6.1 DNA preparation method)
- Pipettes
- Sterile pipette tip
- Vortex mixer
- Centrifuge for microcentrifuge tubes
- Thermal cycler
- Electrophoresis kit
- UV transilluminator

6. PROCEDURE

Please read through the entire procedure before starting.

6.1 DNA Preparation

Various manufacturers offer DNA isolation kits. Please carry out the DNA isolation according to the manufacturer's instructions. The following standard extraction kit is recommended.

Product	Catalog No.	Manufacturer
Viral gene-spin™ Viral DNA/RNA Extraction Kit	17151	iNtRON Biotechnology, Inc.

6.2 Amplification

- ① Prepare appropriate PCR premix tubes and label. And one PCR premix tube for positive control.
- ② Add 2 μ l of template DNA into the PCR premix tube.
- ③ Add 18 μ l of DNase/RNase-free water into the PCR premix tube to total volume as 20 μ l.
- ④ Add 2 μ l of positive control and 18 μ l of RNase-free water into a PCR premix tube for monitoring of amplification and easy interpretation.
- ⑤ Dissolve the blue pellet by pipetting.
Note :The pellet is easily dissolved, by letting the mixture stand at R.T. for 1-2minutes after adding water.
- ⑥ (Optional) Add mineral oil. This step is unnecessary when using a thermal cycler that employs a top heating method (general methods).
- ⑦ Perform PCR reaction of samples as the below process using PCR machine.

PCR cycle		Temp.	Time
1 Cycle	Initial Denaturation	94 °C	5 min.
	Denaturation	94 °C	30 sec.
40 Cycles	Annealing	52 °C	30 sec.
	Extension	72 °C	40 sec.
1 Cycle	Final extension	72 °C	5 min.

6.3 Detection of Amplified Products

- ① Prepare 1.5% agarose gel containing RedSafe™ Nucleic Acid Staining Solution (20,000x).
- ② Load 7 μ l of PCR product and positive control on agarose gel without adding a loading-dye buffer and perform electrophoresis.
- ③ Run electrophoresis by 100V (required about 30~40 minutes).
- ④ Identify the result on ultra-violet (UV) transilluminator.

6.4 Interpretation

- Expected PCR product size : **484 bp**




Fig 1. Electrophoresis of PCR product by **Chicken anemia virus Detection Kit**
 Lane M : 100bp Molecular ladder (iNtRON Biotechnology)
 Lane 1~7 : CAV positive sample
 Lane N : Negative control

6.5 Elimination of carry-over contamination

- Each PCR/RT-PCR Pre-mixture contains 8-methoxypsoralen (8-MOP) for preventing of carry-over contamination.
- All PCR products should be discarded after UV irradiation (10 min/365nm) for preventing from carry-over contamination.

■ NOTE :

- This product was designed to detect more than 100 copies of target gene(or gene segment). When the copy number of target present in the test reaction is less than 100, a false-negative(a negative test result when the attribute for which the subject is being tested actually exists in that subject) may occur. Use this product For Research Use Only.
- Do not use any reagent after the expiration date.
- Do not use together with reagents of other products.
- Follow the instructions.
- Take care in handling of specimen to minimize risk of infection.
- The PCR process is covered by patents issued and applicable in certain countries. iNtRON Biotechnology, Inc. does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is recommended for persons that either have a license to perform PCR or are not required to obtain a license.

- ① In the case of difficult to interpret results due to non-specific bands.
☞ Reduce amount of template by 1/10 dilution and reacts again.
- ② Preparation of PCR reaction at room temperature may cause the non-specific band.
- ③ All procedure should be carried out on ice.

Product	Catalog No.
Chicken Anemia Virus Detection Kit	IP11072
Viral Gene-spin™ Viral DNA/RNA Extraction kit	17151
SiZer™ 100 DNA Marker	24073