

# BCIP/NBT-plus

- Product code:** 3650-10
- Contents:** Contains BCIP 0.692 mM and NBT 0.734 mM in a 2-amino-2-methyl-1-propanol buffer at pH 9.8. Also contains non-toxic stabilizers.
- Applications:** Substrate for alkaline phosphatase conjugates in ELISpot assay.
- Instructions:** Filter the ready-to-use substrate solution (BCIP/NBT-plus) through a 0.45 µm filter and add to ELISpot plates after incubation with alkaline phosphatase and subsequent washing. Develop until spots emerge. Stop the colour development by washing extensively in tap water.
- Storage:** Shipped at ambient temperature. On arrival store substrate refrigerated at 4-8°C. Warm substrate to reaction temperature before use.

**MABTECH**

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Distribuito in ITALIA da  
**Li StarFish S.r.l.**  
Via Cavour, 35  
20063 Cernusco S/N (MI)  
telefono 02-92150794  
fax 02-92157285  
info@listarfish.it  
www.listarfish.it