



DIAGNOSIS

OF RESPIRATORY
INFECTIONS WITHIN
15 MINUTES



HIGH QUALITY RAPID
IMMUNOCHROMATOGRAPHIC TESTS FOR
DIRECT DETECTION OF THE INFECTIOUS AGENT

Sample: swabs, aspirates, lavages



Distribuito in ITALIA da
Li StarFish S.r.l.
Via Cavour, 35
20063 Cernusco S/N (MI)
telefono 02-92150794
info@listarfish.it
www.listarfish.it

> VIRUS

Rapid-VIDITEST Influenza A

Rapid-VIDITEST Influenza A + B

Rapid-VIDITEST Adenovirus Resp.

Rapid-VIDITEST RSV

Rapid-VIDITEST RSV-Adeno Resp.

> BACTERIA


Rapid-VIDITEST Strep A

Rapid-VIDITEST Legionella

Rapid-VIDITEST Streptococcus pneumoniae


Test procedure

> SWABS




> NASOPHARYNGEAL SAMPLE

Add 6 drops of nasopharyngeal wash or aspirate and 9 drops of diluent



> CASSETTES



Wait 10 minutes


> INTERPRETATION OF RESULTS




Positive



Negative



Invalid



Invalid



EUROPEAN UNION
European Regional Development Fund
Operational Programme Enterprise
and Innovations for Competitiveness

* The test procedure is different for Strep A, Legionella and Streptococcus pneumoniae tests. See the instruction manual.