



Li StarFish S.r.l.
 Via Cavour, 35 - 20063 Cernusco S/N (MI), Italy
 Tel. +39-02-92150794 - Fax. +39-02-92157285
 info@listarfish.it - www.listarfish.it

Sensitivities of the serial 4000 assays

Parameter	Catalogue no.	Tests	Sensitivity
Epinephrine (Adrenaline)	LDN-BA E-4100	96	7 pg/mL
Norepinephrine (Noradrenaline)	LDN-BA E-4200	96	35 pg/mL
Dopamine	LDN-BA E-4300	96	7 pg/mL
Epinephrine/Norepinephrine	LDN-BA E-4500	2x96	Epinephrine: 7 pg/mL Norepinephrine 35 pg/mL
Epinephrine/Norepinephrine/Dopamine	LDN-BA E-4600	3x96	Epinephrine: 7 pg/mL Norepinephrine 35 pg/mL Dopamine 7 pg/mL

Li StarFish distribuisce:

