

# pNPP ELISA Substrate

<b>Product code:</b>	3652-P10
<b>Contents:</b>	Contains pNPP, p-Nitrophenyl Phosphate, in diethanolamine buffer.
<b>Applications:</b>	Ready-to-use substrate for ELISA assays with alkaline phosphatase.
<b>Instructions:</b>	Add the substrate after incubation with alkaline phosphatase and subsequent washing. Incubate for typically 30-90 minutes at room temperature. The substrate produces a yellow product. Read the absorbance at 405-410 nm.
<b>Storage:</b>	Shipped at ambient temperature. On arrival store substrate refrigerated at 4-8°C.
<b>Size:</b>	120 ml



Distribuito in ITALIA da  
**Li StarFish S.r.l.**  
Via Cavour, 35  
20063 Cernusco S/N (MI)  
telefono 02-92150794  
fax 02-92157285  
info@listarfish.it  
www.listarfish.it

Note; for research use only.  
Mabtech shall not be liable for the use or handling  
of the product or for consequential, special,  
indirect or incidental damages therefrom.

