



Bio-X Diagnostics S.A.



Cat Nb.	Product	Pathogen(s)	tests/kit	Species	Additional Info
<b>Clostridium Difficile</b>					
BIO K 400	DipFit Clostridium difficile (10x) Combined kits : BIO K 351, BIO K 360, BIO K 374, BIO K 326/2 etc....	10 rapid test strips	ag det.	all species	
<b>Clostridium perfringens</b>					
BIO K 266/2	Monoscreen AgELISA C. perf. alpha toxin/sandwich, double wells	Alpha Toxin of Clostridium perfringens	2-plate kit, 96 samples (strip)	all species	ag det.
BIO K 291/1	Monoscreen AbELISA C. perf. alpha toxin / blocking	Alpha Toxin of Clostridium perfringens	1-plate kit, 96 samples (strip)	all species	Sero blocking
BIO K 291/2	Monoscreen AbELISA C. perf. alpha toxin / blocking	Alpha Toxin of Clostridium perfringens	2-plate kit, 192 samples (strip)	all species	Sero blocking
BIO K 267/1 *	Monoscreen AgELISA C. perf. beta toxin / sandwich, double wells *	Beta Toxin of Clostridium perfringens	1-plate kit, 48 samples (strip)	all species	ag det.
BIO K 267/2	Monoscreen AgELISA C. perf. beta toxin / sandwich, double wells	Beta toxin of Clostridium perfringens	2-plate kit, 96 samples (strip)	all species	ag det.
BIO K 317/2	Monoscreen AbELISA C. perf. beta toxin / blocking	Beta Toxin of Clostridium perfringens	2-plate kit, 192 samples (strip)	all species	Sero blocking
BIO K 268/2	Monoscreen AgELISA C. perf. epsilon toxin / sandwich, double wells	Epsilon toxin of Clostridium perfringens	2-plate kit, 96 samples (strip)	all species	ag det.
BIO K 222/1	Monoscreen AbELISA C. perf. epsilon toxin / blocking	Epsilon toxin of Clostridium perfringens	1-plate kit, 96 samples (strip)	all species	Sero blocking
BIO K 222/2	Monoscreen AbELISA C. perf. epsilon toxin / blocking	Epsilon toxin of Clostridium perfringens	2-plate kit, 192 samples (strip)	all species	Sero blocking
BIO K 269/2	Monoscreen AgELISA C. perf. / sandwich, double wells List of combined kits : BIO K 270/2, BIO K 305/2, BIO K 326/2	Clostridium perfringens	2-plate kit, 96 samples (strip)	all species	ag det.
BIO K 388	DipFit Clostridium perfringens epsilon toxin	10 rapid test strips	ag det.	all species	
BIO K 389	DipFit Clostridium perfringens	10 rapid test strips	ag det.	all species	
BIO K 402	SmartStrips Clostridium perfringens Combined kits : BIO K 306, BIO K 351, BIO K 360, BIO K 372, BIO K 374, BIO K 393 etc....	10 rapid test strips	ag det.	all species	
<b>COMBINED KITS</b>					
BIO K 270/2	Multiscreen AgELISA Enterotoxemia/sandwich, double wells (V2.0)	Alpha, Beta, Epsilon Toxin and Clostridium perf.	2-plate kit, 24 samples (strip)	all species	ag det.
BIO K 326/2	Multiscreen AgELISA Piglet digestive / sandwich, double wells	Rota, E.coli F4, F5, F18, F41, C.diff, C.perf, Crypto	2-plate kit, 12 samples (strip)	pigs	ag det.
BIO K 382/2	Multiscreen AgELISA Alpha toxin and C. perf.	Alpha Toxin and Clostridium perfringens	2-plate kit, 48 samples (strip)	all species	ag det.
BIO K 306	Rainbow Calf Scours (Rota, Corona, E.coli F5 (K99), Crypto, C. perfringens)	5 devices of 5 strips	ag det.	bovines	
BIO K 316	Rainbow Lamb and Kid Scours (Rota, E.coli F5 (K99), Crypto, C. perfringens, Epsilon toxin)	5 devices of 5 strips	ag det.	small ruminants	
BIO K 360	Rainbow Lamb and Kid Scours (Rota, F5, Crypto, C. perfringens, C. difficile)	5 devices of 5 strips	ag det.	small ruminants	
BIO K 372 *	Rainbow Piglet Scours (Rota, F4, F18, C. perfringens) *	5 devices of 4 strips	ag. det.	pigs	
BIO K 374	Rainbow Piglet Scours (Rota, F4, F5, F18, F41, C. perfringens, C. difficile) *	5 devices of 7 strips	ag. det.	pigs	
BIO K 377 *	Rainbow Calf Scours (Rota, Corona, E.coli F5, F41, CS31A, Crypto, C.perf) *	5 devices of 7 strips	ag. det.	bovines	
BIO K 393 *	Rainbow Piglet Scours (Rota, F4, F5, C. perfringens) *	5 devices of 4 strips	ag. det.	pigs	

**NEW PRODUCTS IN RED** \* means: production on request; at least 5 ELISA kits or 5 vials at the time; delivery delay at least 3 working weeks

Comments : ag. det. = antigen detection sec. Ab = secondary antibodies

ind. sero. = indirect serology n.c. = not to be continued \* on request : minimal quantities